





# Automatic Kjeldahl Nitrogen / Protein Analyser


From Traditional to Transformational – Turning Manual Methods into Smart Systems





Labquest, the distinguished line of laboratory equipment from Borosil Scientific Limited, stands as the benchmark of precision and innovation in India's scientific landscape for over six decades. Recognized for unparalleled quality, accuracy, and dependability, Labquest has emerged as the preferred choice in laboratories nationwide. Specifically designed for critical proximate analysis parameters such as protein, moisture, fat, and fiber, Labquest caters to the exacting standards of the food and feed industries.

Our commitment is evident with our specialized services to premier institutions like the Indian Institutes of Technology (IITs). These esteemed academic establishments, renowned for their rigorous research and academic excellence, are our premium customers. Labquest's cutting-edge technologies and indigenously manufactured equipment align seamlessly with the high standards set by IITs, reinforcing our commitment to supporting their scientific endeavors.



## Proven experts for Nitrogen Estimation



Food/Feed/ Dairy

Nitrogen is present in the form of Protein in food products like dairy, vegetables, fruits, meat, fish. Protein perform basic functions of the body such as growth and maintenance, production of enzymes that help in digestion, energy production, blood clotting, muscle contraction. Some proteins are hormones, which are chemical messengers that aid communication between cells, tissues and organs



Agriculture/Soil/ Fertilizer

Nitrates in the soil are used by plants to prepare food. Nitrogen in air can be converted into useful nitrate compounds by bacteria, algae, and even lightning. Once in the soil, nitrogen becomes biologically accessible. This is a major component of chlorophyll, the compound by which plants use sunlight energy to produce sugars from water and carbon dioxide (i.e., photosynthesis).



Wastewater/Sludge

Large amounts of Nitrogen in wastewater utilizes the oxygen released from oxidation of organic waste to form inorganic nitrates. These inorganic nitrates cause marine death. Hence an adequate level of nitrogen is desired in wastewater and the same is detected during wastewater treatment.

#### Others



**Breweries** 



**Pharmaceuticals** 



Rubber/ Dyes

Nitrogen is an essential element required for yeast growth. Individual amino acids, small peptides and some ammonium ions produced from the proteolysis of proteins in barley malt are the main nitrogen sources for yeast metabolism in wort (the liquid extracted from the mashing process during the brewing of alcoholic beverages).

Peptides, enzymes, hormones are all forms of protein that are important constituents for human growth and development.

Peptides are necessary to achieve homeostasis and maintain health, Enzymes are biological catalysts used to speed up reactions involved in many processes like digestion of food, production of energy and synthesis of cellular components and Hormones are involved in blood glucose regulation, appetite regulation



## Digesters

Safe, Simple & Efficient Digestion systems for every kind of application





# Kjeldahl Oracle Digesters



Reproducible digestion due to refined thermal homogeneity. Printer portal for documentation of sample time, temperature and other process parameters with print and PC data transfer ports. Delayed scrubber start ensures safe unattended operation and complete fume evacuation. Cooling stations increase throughput by decreasing the cooling time post digestion.



## Kjeldahl Magnus Digesters

Automatic Kjeldahl digestion unit with 8 / 12 / 20 heating positions of aluminum heating block offers excellent thermal homogeneity and stability

The digesters are easy to use and perform where the user can prepare the samples and leave them unattended after setting the temperature and time. Sulfuric acid used for digestion helps in breaking the bonds of nitrogen. An adequate concentration of catalysts and salt is added to speed up the digestion process. The scrubber system removes the sulfur dioxide fumes keeping the salt: acid ratio undisturbed. Bumping or foaming of samples in tubes is nature dependent and can be avoided by the addition of glass beads or a few drops of hydrogen peroxide.



250 ml straight sample tube provided with the unit is used for sample preparation



Load the manifold for evacuation of digestion fumes



Select the preferred program from data library



Cooled samples are ready for distillation process



Magnus block digester takes over the documentation process from you.



The tube rack with manifold is automatically lowered into the digester and can be lifted automatically in between to check the completion of digestion

The digestion process involves release of pungent sulfur dioxide fumes. Hence the digester is recommended to the use with a fume scrubber only.



## Magnus & Oracle Digester

## Technical Specifications

Model Name	KBD083	KBD123	KBD203	KBD062	KBD082	KBD122
No. of positions	6	12	20	6	8	12
Tube capacity (ml)	250					
Temperature Range	Ambient to 470°C					
Timer Range	1-300 mins					
Temperature Precision	±1°C					
Ramp and Hold	upto 4					
Presets	upto 50					
Auto scrubber ON/OFF	Yes					
PID Autolift	Yes No					
Auto OFF after achieving set time	Yes					
Temperature Calibration	Yes, single and two point					
PC & Printer Connectivity	Yes					
Display	Graphical LCD					
Body	Stainless steel with GRP control pannel					
Current	6A	8A	12A	6A	6A	8A
Power rating	1400W	1600W	2800W	1400W	1400W	1600W
Dimensions (W × D × H)	427 x 490 x 415 x 560 x 415 x 635 x 740 427 x 490 x 665		415 x 635 x 740			
Ordering Code	BLFAKBD083 00000000	BLFAKBD123 00000000	BLFAKBD203 00000000	BLFAKBD062 00000000	BLFAKBD082 00000000	BLFAKBD122 00000000



### Kjeldahl

## Magnus Distillation System



Steam distillation compatible for Nitrogen estimation with sophisticated features makes distillation easy. Easy step by step user modification helps to determine Nitrogen content in the sample.

Soft controls allow the user to manage addition of alkali, boric acid, distilled water and steam with all the other processes being monitored automatically.

Traditionally done using an Allihn's condenser mounted with a Kjeldahl trap, leads to insufficient cooling and extraction of nitrogen and subsequent protein.

The steam distillation using a Friedrichs condenser provides more surface area for cooling. Steam purged into the sample allows thorough mixing for complete nitrogen extraction.



## Kjeldahl Magnus Distillation System

Smart distillation systems for labs. Excellent value for money with high reliability and advance performance.

 Sophisticated construction protects against involuntary contact with hazards acids and chemicals.



Easy and quick Kjeldahl distillation
 Automated addition of chemicals, distillation
 time, steam output with a soft start





 Safety a top priority with tube roader, boiler reagent and door level sensor



A drip tray that protects unit interior



### Kjeldahl

## Oracle Distillation System

Smart, efficient distillation systems primarily for kjeldahl analysis with simple and reliable features.



\*Accessories in the image may not be part of the stranded equipment



### Kjeldahl

## Oracle Distillation System

Single knob control function.

Single press precise addition of alkali anytime during the distillation cycle.



Detachable bottom drip tray for extra protection against acidic attack during sample loading.





- Tube sensor that enables safe operating conditions in event of tube damage.
- Process kill switch to halt the process in case of accidental door opening.



- Reagent refill level sensor to avert reagent dry out.
- Acid and alkali pump calibration for accurate reagent dosing.
- Documentation of no. of cycles and hours of usage with auto service and maintenance alerts.



## Magnus & Oracle Distillation

## Technical Specifications

Feature/ Model	KDI050	KDI040	
Steam addition	Automatic		
Distilled water addition	Automatic		
Alkali addition	Automatic		
Boric acid addition	Automatic		
Tube Drain	Automatic	NA	
Display	7'' colored touch screen	LCD	
Tube Capacity (ml)	250 / 4	400ml	
Steam power (Adjustable)	25-100%	10-100%	
Minimum distillation time (mins)	5	)	
Reagent level sensor	Υe	es	
Anytime alkali addition	Yes		
Tube error detection	Yes		
Door open sensor	Yes		
Process resume from stop	Yes		
One point alkali and distilled water	Yes		
Two point calibration for acid pump	Yes		
Real time monitoring of process parameter	Yes		
SELV protection	Yes		
Programs	pre set as well as user editable		
PC Connectivity	40 preset, 30 user editable	>99.5%	
Nitrogen recovery	>99.5%		
Nitrogen reproducibility	±1 %		
RSD	<1%		
LOD	1-200mg with of ≡0.1mg of Nitrogen		
Current	10A 5.5A		
Power rating	1500 W 1000 W		
Dimension (L × B × H) mm	(480 × 470 × 820) mm (650 × 350 × 725) mm		
Ordering Code	BLFAKDI0500000000 BLFAKDI0400000000		



## Fume Scrubber

Safe, zero exposure & high capacity fume scrubber for Kjeldahl digestion systems





## Fume Scrubber Four Steps for a lab free of toxic fumes



#### SAFE

- Practically zero exposure of harmful sulfur dioxide fumes to user
- Maximum safety with a selfcontained design
- High capacity pump for oil free vacuum pump for larger capacity digesters




#### **EFFICIENT**

- Full control of the Scrubber via the digester
- Maximal neutralization efficiency by the vacuum pump
- Easy compatibility with all digesters

#### Neutralization steps

- Stage 1: Fume extraction and condensation
- Stage 2: Neutralization by distilled water
- Stage 3: Neutralization by alkali
- Stage 4: Neutralization by carbon filter



Model Name	KSC010	KSC020	
Current Consumption	1A	1A	
Pump delivery Rate	28 Ltr/Min	45 Ltr/Min	
Dimensions (WxDxH)	365 X 375 X 480 mm	365 X 375 X 480 mm	
Ordering Code	BLFAKSC01000000000	BLFAKSC0200000000	
HSN Code	84141000		

#### **ACCESSORIES**

Description	Product Code
KSC010 Condenser V2 ASM	BLMAFUTV2KSCCD2932
KSC010 Carbon Bottle ASM	BLM0BGW25012200211
BL 2L Clear Bottle 34/35 Ground Neck ASM	BLA0BGW10001500030
Viton Rubber Tube - 10 mm ID x 14 mm OD	BLA0VRT10D0014DBLK



The sulfur dioxide fumes, a by-product of digestion, causes environmental pollution and is a health hazard if disposed of ineffectively. The need for a fume hood can be nullified by the fume scrubber system by LabQuest Borosil.



The easy to assemble construction makes all components easily accessible

The protective front shield gives extra protection to the scrubber body



The scrubber helps neutralizing toxic fumes by its 4 stage mechanism

Stage 1 : Condensation via condenser Stage 2 : Neutralization by distilled water

Stage 3: Neutralization by Alkali

Stage 4: Neutralization by carbon filter

#### Technical Data & Ordering Info

Model Name	KSC010	KSC020
Current	1A	1A
Pump delivery Rate	28 Ltr/Min	45 Ltr/Min
Dimensions (WxDxH)	(445 × 375 × 460) mm	
Ordering Code	BLFAKSC0100000000	BLFAKSC0200000000

#### Accessories

Description	Product Code
Silicon Transt Tube-8MMIDX10MMOD - 1Mtr	BLA0STT8D00012TRNP
Viton Rubber Tube-10MMIDx14MMOD - 1Mtr	BLAOVRT10D0014DBLK
KSC Condenser V2 ASM	BLAAFUTV2KSCCD2932
KSC Carbon Bottle ASM	BLAACRB02020018000
2L Reagent bottle - 1500030 ASM	BLA0BGW10001500030



### **Application Center**

Navigating the ever-evolving terrain of scientific exploration, our Application Development Center, strategically located at our manufacturing hub in Maharashtra, serves as a dedicated resource ready to assist you. Our center offers comprehensive support in method validation, method development, sample testing, and all other aspects of application or technical needs.

Our dedicated team at the Application Development Center understands the unique challenges faced by IITs and stands ready to collaborate with you. Whether it's refining existing methodologies or developing new ones, we aim to enhance your speed, efficiency, and productivity in the laboratory. By providing tailored solutions, we ensure that your evolving needs are met, setting the stage for optimal performance and success.

Partner with us at the forefront of innovation, where the collaborative spirit of our Application Development Center meets the pioneering efforts of IITs. Together, let's navigate the frontiers of scientific exploration and achieve new heights of excellence in research and discovery.





## Do's and Don'ts

## For Nitrogen estimation system

Do's	Don'ts	
	Digester	
Unit should be operated by trained personnel only	Do not touch the hot top surface while it is working and immediately after completion of digestion	
It is ideal to use an adequately weighed and pre mixed mixture of catalyst and salts	No blank tubes to be used. If fewer samples are to be processed, other tubes are to be operated with blank	
Check the power supply with proper earthing	Do not operate the unit beyond below mention load/day 6 position: 18 samples/day 8 position: 24 samples/day 12 position: 36 samples/day 20 position: 60 samples/day	
Clean the suction end of every day after use	While process do not remove fume extractor tray	
Clean the suction hose from manifold to condenser end of every day after use		
Keep all test tube with test tube tray		
Making sure the tabs provided in the fume extractor tray gets fitted into the slot provided on the test tube tray.		
	Scrubber	
Distilled water should replace the solution after every digestion cycle.	Do not operate the unit without NaOH solution or fully neutralized solution. This may damage the pump diaphragm.	
While process open the condenser inlet water	End of the day close the tap water	
NaOH in the scrubber should be changed once it turns muddy black or after checking the acidity with phenolphthalein indicator	Do not stop scrubbing for at least 15 mins post completion of digestion time	
Distillation unit		
Unit should be operated by trained people only	Do not open the door while in the distillation process	
Press the knob / touch the screen gently while using the unit	Do not try to open the back door of the unit	
Check the power supply with proper earthing	Do not interchange the reagent tubes	
Always check the distilled water TDS before preparing chemicals	Do not operate the unit beyond 60 samples in a 8 hrs shift day	
Dip sensors properly in the reagent bottles	Do not touch the hot surface wherever written on the unit	
Gently open the door of the unit	Do not try to remove control knob from the unit	
Always check the water supply and pressure of the water tap		
All the inlet pipes should be connected properly		
Always clean the NaoH pipe & instrument after analysis		
Priming the pumps and hoses with the distilled water everyday before commencing analysis		
Exchange acids and alkalis for distilled water when not in use		



## Service capabilities

At Borosil Scientific Limited, we understand that for premier institutions like the Indian Institutes of Technology (IITs), seamless laboratory operations are paramount.

Our proactive approach ensures that your Labquest instrument enjoys an extended lifespan with minimal maintenance and service expenses throughout its entire lifecycle. A dedicated team of engineers stationed in Kolkata and Guwahati stands ready to provide swift and seamless support, addressing any issues promptly.

To further enhance your experience, our strategically stocked spare inventory at the regional warehouse in Kolkata guarantees part availability in less than 48 hours. This ensures that your laboratory operations at IITs remain uninterrupted, minimizing downtime and optimizing efficiency.

With Borosil Scientific, your Labquest instrument transcends being just a purchase; it becomes a strategic investment backed by our unwavering commitment to lasting performance and unparalleled service.







#### Digester Accessories

#### Cooling stations for faster sample cooling and increased sample throughput



**BLGOFUTKBD** 080D2510 Fume Extractor



**BLGAFUT** 100PTFE1651 Teflon Cap



**BLGATTBO** 1041250001 Digestion Tube



6 POSITION



20 POSITION

#### Description

Viton Rubber Tube -10MMIDx14MMOD - 1Mtr BLA0VRT10D0014DBLK Teflon Fume Cap Accs (Case Lot -01 No) KBD062 Test Tube Holder Accs KBD082 Test Tube Holder Accs KBD122 Test Tube Holder Accs KBD062 Fume Extractor Accs KBD082 Fume Extractor Accs KBD122 Fume Extractor Accs KBD060 Fume Extractor Glass Accs KBD080 Fume Extractor Glass Accs KBD062 Cooling Tray Asm Accs KBD082 Cooling Tray Asm Accs

#### **Product Code**

BLAAFUT100PTFE1651 BLAAKBD06220000000 BLAAKBD08220000000 BLAAKBD12220000000 BLAAKBD06230000000 BLAAKBD08230000000 BLAAKBD12230000000 BLA0FUTKBD060D2510 BLA0FUTKBD080D2510 BLAAKBD06170000000 BLAAKBD08250000000

#### Description

KBD123 Cooling Tray Asm Accs KBD203 Cooling Tray Asm Accs KBD083 Test Tube Holder Accs (Autolift) KBD123 Test Tube Holder Accs (Autolift) KBD203 Test Tube Holder Accs (Autolift) KBD083 Fume Extractor Accs (Autolift) KBD123 Fume Extractor Accs (Autolift) KBD 203 Fume Extractor Accs (Autolift) KBD200 Fume Extractor Glass Accs 250 ml Borosilicate Test Tube ASM (Case lot-5) BLMATTB01041250000 Quartz Test Tube 250ml ASM (Case lot-2) 400 ml Borosilicate Test Tube ASM (Case lot-5) BLAATTB01050400001

#### Product Code

BLAAKBD12350000000 BLAAKBD20350000000 BLAAKBD08320000000 BLAAKBD12320000000 BLAAKBD20320000000 BLAAKBD08330000000 BLAAKBD12330000000 BLAAKBD20330000000 BLA0FUTKBD200D2510 BLGATTQ01041250000

#### Combine digester with a Bottle top dispenser for safe acid dispensing

Bottle Top Dispensers offer simple solution for safe dispensing of reagents & chemicals in desired volume.

#### Distillation Acessories



Condensor

Dosing head

#### Description

Silicon Transt Tube - 8MMIDX10MMOD - 1Mtr Clear Braided PVC Hose Pipe ACCS - 1Mtr KDI PTFE Tube - 10mm OD X 08mm ID - 1Mtr Dosing head Condensor KDI Flow Switch ASM ACCS **KDI Drip Tray ASM** 

Oven Mittes - KDI ASM KDI PP Jerry Can ACCS

#### Product Code

BLA0STT8D00012TRNP BLA0PIPPVC0HOS00CB BLA0TUB000PTFE1008 BLG0DOH050KJV2DIST BLG0CON050KJV2DIST BLAAKDI05010000001 **BLAOKDIOOODRPOTRAY** BLA0OVM000KDI00000 BLA0JCNPP00000005L



## Care and maintenance

### Of Nitrogen estimation system

Follow the change time interval of consumables as given below (to be changed by service engineer only)

	Component	Change Time Interval
1	Silicon Tube	6 months
2	Silicon Tube for Peristaltic Pump	3-6 Months condition based.
3	Viton Tube	3-6 Months condition based.
4	Braided Hose	Yearly
5	Teflon Pipe	Yearly
6	GL connectors	Yearly
7	Glass parts ( Dosing Head /Condenser )	Yearly
8	NRV	Condition based
9	Solenoid valves	Condition based
10	Heaters	Condition based
11	Drain Valves	Condition based
12	Limit / Flow switches	Condition based
13	Level Sensors	Condition based
14	PT100 Temperature Sensor	Condition based
15	Dosing head holder	6 Months condition based.
16	Pumps	Condition based



1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051. India T :+91 22 6740 6300 F :+91 22 6740 6514 E : borosil@borosil.com W : www.borosil.com

Pune Manufacturing Unit: BOROSIL SCIENTIFIC LTD. 21-24, Plot No. 8, Indialand Global Industrial Park, Hinjewadi Phase - 1, Pune - 411 057







✓ lab.support @borosil.com **\** 1800-224-551