



# Hybrid e- Dispenser

The Ultimate Laboratory Dispensing solution

# Top-Notch Features for Precision and Accuracy

Precision Control  
Knob for effortlessly  
toggling between 3  
preset speeds.


Distinct Modes  
for Dispensing  
and Titration.

Accuracy and  
Precision Complying  
with ISO 8655

USB interface  
facilitates transfer  
of data to PC

QR scanner for easy  
access to Chemical  
compatibility chart

Includes a range of  
adapters for  
mounting on diverse  
bottle neck



SMART

RUGGED

PRACTICAL



Reliable results with no  
downtime with inbuilt  
calibration feature



Select between the 3 preset  
speed- rapid dispensing ,  
intermediate dispensing and  
meticulous dispensing for  
precise control



Amber colored and transparent  
Inspection window



Enhanced Visibility - Large and  
clear display ensuring visible  
reading from every angle.

# Technical Specifications

Volume	50mL
Volume Accuracy in Dispense mode	+/- 0.30ml
Volume Accuracy in Titration mode	+/- 0.05ml
Dispense Incremental Value	0.2mL
Dispenser Modes	Time based, Fixed volume dispensing
Speed for Titration	3 Calibrated Pre-set Speeds
Operation Mode	Dispensing + Titration
Viewing Window	Amber & Clear Transparent
Operating method	Piston Operated
Inbuilt Calibration	Yes
Adapters	GL28, GL32, GL36, GL38, GL40 & GL45
Connectivity for data Storage	USB interface for PC connectivity, enables storage of upto 10 readings



**BOROSIL®**  
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,  
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India.  
T +91 22 6740 6300 F +91 22 6740 6514  
E bsl@borosil.com W www.borosilscientific.com

PROUDLY MADE IN   
india

