



Rapid Fat Analyser

**Distinctive Fat Estimation Systems
by Randall Extraction**

Technical Specifications

Model Name	FAS060
No of positions	6 No's
Sample size	0.5 gm to 10 gm
Recovery of Solvent	80-90%
Limit of detection	0.1% to 100 %
Heating profile	Energy regulator 1 to 10
Display	7" touch screen
Water flow sensor	Yes
Beaker solvent capacity	150 ml
Condenser	Seal Viton (FKM)
Chiller water temp	10-15 degree
Heater temp range	280 degree Max
Recovery bottle capacity	2 liter
Flow cut off sensor	yes
Vacuum pump flow rate	380 L/hr
Power supply	230 V 50 – 60 Hz
Current Consumption	5.2A
Power consumption	1200 W
Dimensions (LxBxH) mm	840 x 430 x 600
HS Code	84198990
Ordering code	BLFRFAS060000000000
Price/ Piece	On request

BOROSIL®
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India.
T +91 22 6740 6300 F +91 22 6740 6514
E bsl@borosil.com W www.borosilscientific.com

PROUDLY MADE IN
india

