



# HYBRID E - DISPENSER

The Ultimate Laboratory Dispensing solution

# Top-Notch Features for Precision and Accuracy

Precision Control  
Knob for effortlessly  
toggling between 3  
preset speeds.


Distinct Modes  
for Dispensing  
and Titration.

Accuracy and  
Precision Complying  
with ISO 8655

USB interface  
facilitates transfer  
of data to PC

QR scanner for easy  
access to Chemical  
compatibility chart

Includes a range of  
adapters for  
mounting on diverse  
bottle neck





Reliable results with no  
downtime with inbuilt  
calibration feature



Select between the 3 preset  
speed- rapid dispensing ,  
intermediate dispensing and  
meticulous dispensing for  
precise control



Amber colored and transparent  
Inspection window



Enhanced Visibility - Large and  
clear display ensuring visible  
reading from every angle.

# Technical Specifications

|                                   |  |
|-----------------------------------|--|
| Volume                            | 50mL   |
| Volume Accuracy in Dispense mode  | +/- 0.30ml   |
| Volume Accuracy in Titration mode | +/- 0.05ml   |
| Dispense Incremental Value        | 0.2mL  |
| Dispenser Modes                   | Time based, Fixed volume dispensing                                    |
| Speed for Titration               | 3 Calibrated Pre-set Speeds  |
| Operation Mode                    | Dispensing + Titration   |
| Viewing Window                    | Amber & Clear Transparent  |
| Operating method                  | Piston Operated  |
| Inbuilt Calibration               | Yes  |
| Adapters                          | GL28, GL32, GL36, GL38, GL40 & GL45                                    |
| Connectivity for data Storage     | USB interface for PC connectivity, enables storage of upto 10 readings |



Visit Website



Connect with us  
on LinkedIn



**BOROSIL®**  
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,  
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India  
T : +91 22 6740 6300 F : +91 22 6740 6514  
E : borosil@borosil.com W : www.borosil.com

Pune Manufacturing Unit:  
**BOROSIL SCIENTIFIC LTD.**  
21-24, Plot No. 8, Indialand Global Industrial Park,  
Hinjewadi Phase - 1, Pune - 411 057

**Service Support:**

✉ lab.support@borosil.com  
☎ 1800-224-551