



MAGNETIC STIRRER

One Compact Solution for Reliable Stirring

PRODUCT SPECIFICATION

PARAMETERS	MGS001
Stirring Position	1
Max Capacity	2L
Speed Range	300 - 2200 RPM
Speed Control	PID
Speed Feedback	Yes
Display	2 X LED 7 Segment
Controlling	Rotary Knobs & Push Keys
Timer	0 - 999 min & Infinite
Top Plate Material	SS304
Top Plate Dimension (LxB)	155mm X 145mm
Product Dimensions (LxBxH)	285mm x 190mm x 120mm
Magnetic Stir Bar	25mm
Fuse Rating	0.5 Amps
Power Consumption	0.5A, 230V, 25W



SAFE



RELIABLE



RUGGED

Visit Website



Connect with us
on LinkedIn



BOROSIL®
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Pune Manufacturing Unit:
BOROSIL SCIENTIFIC LTD.
21-24, Plot No. 8, Indialand Global Industrial Park,
Hinjewadi Phase - 1, Pune - 411 057

Service Support:

✉ lab.support@borosil.com
☎ 1800-224-551