




DRY BATH

Dynamic | Safe | Simple


Tapti (ताप्ती)
Dry Bath




- LabQuest digital dry bath with aluminum blocks provides a rapid and controlled heating environment in micro centrifuge tubes.
- The high-wattage heater is microprocessor controlled for accurate temperature control.
- Temperature is set by adjusting the easy-to-operate rotary knob.
- Interchangeable blocks for requirement of every lab.

Model	LBH010
Temperature Accuracy	±0.5°C at 37°C
Display Resolution	±0.1°C
Heating Up Time	≤ 18 min 30°C TO 130°C
Process Timer	0-999 min
Temperature Control	PID
External Probe	PT 100 included
Internal Probe	PT 100
Unit External Dimensions (W x D x H)	(285 x 190 x 156) mm
Cover for Dry Block	UV Protection
Power Consumption	200 W
Current Consumption	1 A
Aluminum Block	To be order separately
HSN Code	84198990
Ordering Information	BLFHLBH01000000000


- 1 Versatile multi well block** for 8 wells of 0.2ml tube, 8 wells of 0.5ml tube, and 8 wells of 2.0ml or 1.5ml
BLAALBHVMW0008W004



- 2 8 well PCR strip** of 0.2 ml, 32 wells, $\phi 6.35\text{mm}$, depth 19mm
BLAALBHPCR3208W004



- 3 For 0.5 ml tube,**
24 wells, $\phi 8.0\text{mm}$, depth 25mm
BLAALBH0MW2406W004




- 4 FOR 1.5ml or 2.0 ml tube**
16 wells of , $\phi 11.0\text{mm}$, depth 30mm
BLAALBH0MW1604W004




- 5 15 ml tube***
12 wells $\phi 17.3\text{mm}$, depth 50mm
BLAALBH0MW1204W004




- 6 For 50 ml tube***
4 wells, $\phi 29.2\text{mm}$, depth 50mm
BLAALBH0MW0402W004



- 7 Flat 96 well plate block***
BLAALBH0MW0096W004



- 8 Bead bath or water bath block**
Vol of 300ml#
BLAALBH0MW000BB004



#Supplied without Beads

*During usage of 15ml, 50 ml & 96 Well Plate Aluminum Block, Cover remains open

Visit Website



Connect with us
on LinkedIn



Service Support:

✉ lab.support@borosil.com
☎ 1800-224-551