



Glassdryer

Engineered for speed and simplicity, the Labquest Glass Dryer delivers quick, efficient drying of lab glassware in under 15 minutes. Ideal for Kjeldahl tubes, and daily labware, it's a space-saving solution built for convenience—whether on a bench or wall-mounted.

Key Features



Interchangeable Nozzles

Adaptable for all varied glassware types with easy-to-swap drying nozzles (16, 10, 7 mm)



Quick Drying with Auto Shut-Off

Dries glassware at 60 °C in 10–15 mins with set-and-forget timer functionality.



Climate-Adaptive Drying

Customize time and heat for seasonal variations, perfect drying all year round.



Hassle-Free Maintenance

Removable dust filters ensure clean airflow and easy upkeep. It also provides drip collection tray.



Flexible Design

Rugged GI powder-coated body and corrosion-free design for long-term reliability. Can used as benchtop or also as wallmounted

Technical Data & Ordering Information

Model Name	GDS020
Number of Nozzles	21
16 mm	8 Nos
10 mm	8 Nos
7 mm	5 Nos
Temperature Range	Ambient to 70°C
Time Control	Yes
Drying Power	650 W
Dimension (L × B × H) in mm	440 × 170 × 550
Voltage (V)	240 V
HSN Code	84798200
Ordering Code	BLFGGDS02000000000



BOROSIL®
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Pune Manufacturing Unit:
BOROSIL SCIENTIFIC LTD.
21-24, Plot No. 8, Indialand Global Industrial Park,
Hinjewadi Phase - 1, Pune - 411 057

Visit Website



Connect with us
on LinkedIn



Service Support:

✉ lab.support@borosil.com
☎ 1800-224-551