



SOXHLET APPARATUS

Your Distinctive Partner in Fat Estimation Systems

The Only Soxhlet Apparatus With ACTUAL TIME SETTING AND ENERGY PERCENTAGE DISPLAY

Carry out soxhlet extraction up to 6 samples in parallel with the series Gen 2 heating mantles. The functional, robust construction keeps your running costs low and at the same time minimizes the risk of failure. The compact design ensures work safety and ergonomic working

SAFE

Inbuilt safety control



INTELLIGENT

LED display of time, safe temperature and energy percentage



PRECISE


Inbuilt Precise temperature control upto $\pm 1^{\circ}\text{C}$



START YOUR SOXHLET AND WALK AWAY!



The smart SOX series of soxhlet extractors



Technical Sepcifications



SOX012
1 position soxhlet



SOX062
6 position soxhlet

Model Name	SOX012	SOX015	SOX042 *	SOX045 *	SOX062	SOX065
Position	1	1	4	4	6	6
Extractor capacity (ml)	100	200	100	200	100	200
Flask Capacity (ml)	250	500	250	500	250	500
Temperature Accuracy	±1°C					
Mantle Temperature	Upto to 320°C					
Safe temperature setting	Yes					
Knob for Extractor	Yes					
Display Timer	1-999 min & Continuous					
Dimension (W x D x H)mm	200 x 300 x 800		750 x 300 x 800		1050 x 300 x 800	
Weight (Kgs)	7		12		15	
Power (W)	400		1600		2400	
HSN Code	84194090					
Ordering code	BLFRSOX0120 0000000	BLFRSOX0150 0000000	BLFRSOX0420 0000000	BLFRSOX0450 0000000	BLFRSOX0620 0000000	BLFRSOX0650 0000000

*Till stock lasts

BOROSIL®
Scientific

1101, Crescenzo, G-Block, Opp. MCA Club,
Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India
T : +91 22 6740 6300 F : +91 22 6740 6514
E : borosil@borosil.com W : www.borosil.com

Pune Manufacturing Unit:
BOROSIL SCIENTIFIC LTD.
21-24, Plot No. 8, Indialand Global Industrial Park,
Hinjewadi Phase - 1, Pune - 411 057

Service Support:

✉ lab.support@borosil.com
☎ 1800-224-551

Visit Website



Connect with us
on LinkedIn

