







Microscope slides & Cover glasses


Reliable Clarity for Every Specimen

Our microscope slides are engineered for clarity & consistency, our microscope slides are ideal for precise sample examination. Crafted from high-quality glass, they ensure smooth edges, excellent transparency, and compatibility with all standard microscopes.

Uses of

Microscope Slides


Microscope slides are essential tools for preparing & examining specimens under a microscope. They are widely used in biological, chemical, and medical laboratories to hold tissue samples, cells, microorganisms, or chemical substances for observation & analysis.



Material

Made of Sodalime glass with thickness of 1.0 - 1.2mm

Types



Packaging

Packed in a paper box containing silica desiccant to absorb moisture 50 slides per packet & 25 packets per case

Product Code	Description	Dimensions (mm) 18 X 18	Quantity Per Packet (No)	Packets / Case
*9100P02	Plain, Ground Edges	76 X 26 X 1	50	25
9100D01	Double Frosted Ground Edges	76 X 26 X 1	50	25


^{*}ISI marking is applicable only for plain packaging with grounded edges.



Scan for **Product Details**

Why Borosil Cover Slips?

- High Optical Clarity Ideal for detailed microscopic observations.
- Uniform Thickness Ensures consistent focus and imaging.
- Cleanroom Packing Reduces contamination risks.
- Shatter Resistance Carefully annealed to reduce breakage during handling.



Square - 9115501



Rectangular - 9115R01



Square - 9115S02



Rectangular - 9115R02



Square - 9115S03



Rectangular - 9115R03



Conformance to Standard ASTM E211

Product Code	Description	Dimensions (mm)	Quantity / Packet	Packets / Case
9115801		18 X 18		
9115802	Square	22 X 22	100	50
9115\$03		24 X 24		
9115R01	Rectangular	22 X 40	100	50
9115R02		22 X 50		
9115R03		22 X 60		



Scan for **Product Details**



CIN: U74999MH1991PLC061851

Visit **Website**



Follow us on **LinkedIn**



Registered & Corporate Office:

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India.

T +91 22 6740 6300 • F +91 22 6740 6514

Write to us for export enquiries at bsl.exports@borosil.com

MUMBAI

mumbaisales@borosil.com

+91 22 6740 6400 +91 22 6740 6444 **DELHI**

delhi@borosil.com +91 11 4150 5893 CHENNAI

chennai@borosil.com +91 44 2822 6012 / 13 +91 44 2822 6014 KOLKATA

calcutta@borosil.com +91 33 2229 9166 All rights reserved. Borosil Scientific Limited. 2025.

Any unauthorized use, copying, or distribution of this brochure is strictly prohibited and constitutes a legal violation. Borosil Scientific Limited reserves all rights to take legal action against individuals or entities who violate these legal rights.