



Certified sampling Storage vials

Ideal for chromatography,
Storage

Made from 5.1 Low Expansion
Borosilicate Glass



Our sampling vials are designed for reliable sample containment, storage, and transport in a variety of analytical, pharmaceutical, & laboratory applications. Crafted with high quality materials, they ensure sample integrity and minimize contamination risk.

How Borosil Storage Vials different from others?

Provide with Quality Certificate

Packed in ISO "Class 8" clean room facility

Complies glass composition & Spectral transmission test as per USP chapter 660/USP 40 & European pharmacopoeia 9.0



Direct standard



Avoiding wastage



Protected Contamination



Air-tight protected



| Capacity ml | Approx Mouth ID mm | Dimensions O. D X Length mm | Case Lot (Nos) |
|-------------|--------------------|-----------------------------|----------------|
| 3 | 10 | 16 X 36 | 1200 |
| 5 | 10 | 16 X 55 | 900 |
| 10 | 16.5 | 24 X 50 | 1050 |
| 15 | 16.5 | 24 X 58 | 900 |
| 20 | 16.5 | 24 X 70 | 750 |
| 30 | 16.5 | 24 X 95 | 600 |

BOROSIL[®]
Scientific

Borosil Scientific Limited

CIN : U74999MH1991PLC061851

Visit Website



Follow us on LinkedIn



Scan for Product Details



Registered & Corporate Office :

1101, Crescenzo, G-Block, Opp. MCA Club, Bandra Kurla Complex, Bandra (E), Mumbai - 400 051, India.

T +91 22 6740 6300 ● F +91 22 6740 6514

E bsl@borosil.com ● W www.borosilscientific.com

Write to us for export enquiries at bsl.exports@borosil.com

MUMBAI

mumbaisales@borosil.com

+91 22 6740 6400
+91 22 6740 6444

DELHI

delhi@borosil.com

+91 11 4150 5893

CHENNAI

chennai@borosil.com

+91 44 2822 6012 / 13
+91 44 2822 6014

KOLKATA

calcutta@borosil.com

+91 33 2229 9166

All rights reserved. Borosil Scientific Limited. 2025.

Any unauthorized use, copying, or distribution of this brochure is strictly prohibited and constitutes a legal violation. Borosil Scientific Limited reserves all rights to take legal action against individuals or entities who violate these legal rights.